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### Specifications

- *Sample diameter  $5 \text{ mm} < D < 90 \text{ mm}$*
- *range of radius of curvature  
flat to  $R > 1 \text{ m}$  (concave) and  
flat to  $R > 1 \text{ m}$  (convex)*
- *rel. precision of measurement  
 $0,1\% @ R = 1 \text{ m}$   
 $1,0\% @ R = 60 \text{ m}$*
- *length of the optical bench: 800 mm*
- *moving range: 600 mm*

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
### Company profile

*Dioptic supplies you with patented diffractive reference objectives (DFNLs) for high-precision interferometric measurements of optical surfaces. Furthermore, Dioptic's highly qualified team of physicists and engineers develops solutions for your problems in the fields of optical measurement technology, photometry, fluorescence photometry, spectroscopy, microscopy, and interferometry. Working in these fields for many years now, Dioptic is very experienced in setting up customized measurement and analytical systems.*

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
## Compact setup for testing large radii of curvature

*The best-possible precision for measuring the absolute radius of curvature is provided by an interferometric measurement. This method, however, is limited to typically  $lR < 1$  m by the increasing length of the optical path.*

*The Large Radius Test provides a continuous measurement range for radii of curvature from  $R = 1$  m concave to flat and from  $R = 1$  m convex to flat (see area marked in blue), with interferometric precision, in a compact and convenient setup.*

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## Continuously adjustable test wavefront

*The combination of two customized measurement objectives, a diverger and a converger, generates a continuously adjustable test wavefront. Its radius of curvature may be changed by the distance between the two objectives. Depending on the direction the converger is moved in, we obtain a weakly curved convergent or divergent test wavefront.*

*According to the required precision of the topography measurement, the wavefront error can be reduced to  $\leq \lambda/5$  for the complete measurement range.*

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## comfortable measurement

*As soon as the test wavefront has been adjusted to the test sample by moving the converger (power=0 in the interferogram), a software calculates the radius of curvature from the positions of the objectives and the test sample. Depending on the customer's needs, the positioning can be implemented motorized or manually.*